

Application Work AW US6-0249-062017

Analysis of Chloride, Phosphate, Malate, Sulfite, Tartrate, Sulfate, and Oxalate in Red and White Wine

Branch

Food, stimulants, beverages, flavors

Keywords

IC; 940; Conductivity; Chloride; Phosphate; Sulfate; Sulfite; Malate; Tartrate; Oxalate; Anions; 858; MagIC Net 3.2; Metrosep A Supp 10-100/4.0; Ultrafiltration; Wine

Summary

Product consistency and quality is of utmost importance to winemakers. It is also critical to monitor and evaluate yeast performance and efficiency throughout the fermentation process. This wine analysis can aid vintners with ensuring predictable flavor and aroma characteristics in finished wine by monitoring common indicators of acidity, mouthfeel, and balance. It also evaluates nutrients and other additives which could potentially have deleterious effects on efficiency and production during the fermentation process. This application work shows the use of Metrohm IC to analyze red and white wine for chloride, phosphate, sulfite, sulfate, malate, tartrate, and oxalate. For this work, a Metrohm Professional IC with sequential suppression and conductivity detection was used in combination with the 858 Sample Processor with Ultrafiltration. The in-line Ultrafiltration protects the column and system components from particulates and ensures trouble-free operation. The rapid analysis enables high throughput laboratories to maximize production. Dosino regeneration and STREAM rinsing are used for the MSM.

Reagents

- Sodium Carbonate, CAS 497-19-8: Sigma Aldrich 223484-500G, 99.95 – 100.05%
- Sodium Bicarbonate, CAS 144-55-8: Sigma Aldrich S6014-500G, 99.7 – 100.3%
- Sodium Sulfite, CAS 7757-83-7: Sigma Aldrich 71989-250G, 98%
- Tartaric Acid, CAS 87-69-4: J.T. Baker 4104-01
- Malic Acid, CAS 617-48-1: J.T. Baker P494-07
- Oxalic Acid, CAS 144-62-7: Sigma Aldrich 658537-100G, 99.999% trace metals basis
- Perchloric Acid, CAS 7801-90-3: Sigma Aldrich 311421-50mL, 99.999% trace metals basis
- Sulfuric acid, CAS 7664-93-9
- Ultrapure water, resistivity >18 MΩ·cm (25 °C)
- Metrohm Individual Anion Standards, 1000 ppm

- Chloride, Metrohm USA ERA-IC1002
- Phosphate, Metrohm USA ERA-IC
- Sulfate, Metrohm USA ERA-IC1007

Instruments

940 Professional IC Vario ONE/SeS	2.940.1400
IC conductivity detector	2.850.9010
858 Professional Sample Processor	2.858.0020
Sample Rack 148 place + 3 Special Beakers	6.2041.440
MagIC Net 3.2	6.6059.322
Metrosep A Supp 10 - 100/4.0	6.1020.410
800 Dosino	2.800.0010
Dosing Unit 2 mL	6.3032.120
Ultrafiltration Equipment	6.5330.110
20 µL sample loop	6.1825.210
MSM Rotor A	6.2832.000
Filtration Membrane, 0.2 µm	6.2714.020
Pump tubing LFL (yellow/yellow)	6.1826.390
Pump tubing LFL (black/black)	6.1826.340



Solutions

Eluent: 5.0 mmol/L sodium carbonate + 5.0 mmol/L sodium bicarbonate + 5 µmol/L perchloric acid

Regenerant for Suppressor: 500 mmol/L sulfuric acid

Samples

Red wine
White wine

Standards

1000 ppm stock standards were used for all analytes. Sulfite, tartrate, malate, and oxalate stocks were prepared

from salts. Chloride, phosphate, and sulfate were prepared using commercially purchased standards.

Sulfite was calibrated separately due to sulfate contamination in the sodium sulfite salt. Sulfite stock standard was fixed with 10% isopropanol and working standards were fixed with 2% isopropanol in order to prevent oxidation of sulfite to sulfate.

In ultrapure water (ppm)

	S1	S2	S3	S4	S5	S6
Chloride	1	5	10	25	50	100
Phosphate	1	5	10	25	50	100
Malate	1	5	10	25	50	100
Sulfite	1	5	10	25	50	100
Tartrate	1	5	10	25	50	100
Sulfate	1	5	10	25	50	100
Oxalate	1	5	10	25	50	100

Sample Preparation

Wine samples were gravimetrically diluted 1:10 and 1:50 in ultrapure water. Vial caps were used to minimize oxidation. Samples were then directly injected by the 858 with ultrafiltration.

IC Parameters

Anion eluent flow	1.0 mL/min
Column temperature	35°C
Injection volume	20 µL
Anion P _{max}	20 MPa
MSM Rinsing	STREAM
MCS	On
Degasser	On
Recording Time	20 minutes

Calculation

Automatic integration with MagIC Net 3.2 software using peak area for all analytes.

Results

All results summarized in the appendix.

Comments

This application may be deployed utilizing a 930 Compact IC Flex, and is also suitable for use with Metrohm's In-Vial Dilution Technique (MiVDT).

Acetate was not quantified and is shown in the sample chromatograms for identification purposes. If oxalate is not present, the run time may be slightly reduced.

References

AW US6-210-122014: Metrohm In-Vial Dilution Technique (MiVDT) for Anions with Intelligent Dilution, Ultrafiltration, Analyte Delimiter Logic, and Eluent Production Module


Date

June 14, 2017

Author


Ashley R. Wittig
Metrohm USA Inc.

Instrumentation Set-Up




Calibration Curves

Chloride (Anions)




Function: $A = -0.0264217 + 0.0118155 \times Q + 3.20010E-8 \times Q^2$
 Relative standard deviation 0.261555 %
 Correlation coefficient 0.999999

Phosphate (Anions)




Function: $A = -8.89112E-3 + 3.66144E-3 \times Q + 6.57402E-8 \times Q^2$
 Relative standard deviation 0.843057 %
 Correlation coefficient 0.999985

Malate (Anions)




Function: $A = -7.98903E-4 + 3.30240E-3 \times Q - 2.78528E-7 \times Q^2$
 Relative standard deviation 1.603212 %
 Correlation coefficient 0.999939

Sulfite (Anions)




Function: $A = -0.0340779 + 3.81770E-3 \times Q + 2.23093E-7 \times Q^2$
 Relative standard deviation 1.366768 %
 Correlation coefficient 0.999964

Tartrate (Anions)




Function: $A = -1.87955E-3 + 4.33315E-3 \times Q - 2.67851E-7 \times Q^2$
 Relative standard deviation 0.476093 %
 Correlation coefficient 0.999995

Sulfate (Anions)



Function: $A = -0.0155494 + 8.35366E-3 \times Q + 1.34964E-7 \times Q^2$
 Relative standard deviation 0.621130 %
 Correlation coefficient 0.999992

Oxalate (Anions)



Function: $A = -3.78364E-3 + 7.29359E-3 \times Q - 3.80156E-7 \times Q^2$
 Relative standard deviation 0.884738 %
 Correlation coefficient 0.999982


Results

Red Table Wine (ppm)	Chloride	Phosphate	Malate	Sulfite	Tartrate	Sulfate	Oxalate	Dilution
Replicate 1	60.0	771.9	92.2	26.1	1754.5	553.0	<10	10 (50 for tartrate)
Replicate 2	60.0	769.7	92.1	26.9	1754.9	553.0	<10	10 (50 for tartrate)
Replicate 3	60.0	773.3	92.1	27.0	1757.5	553.0	<10	10 (50 for tartrate)
Average	60.0	771.6	92.1	26.7	1755.6	553.0	<10	
RSD (%)	0.027	0.23	0.064	1.9	0.092	0.007		


White Table Wine (ppm)	Chloride	Phosphate	Malate	Sulfite	Tartrate	Sulfate	Oxalate	Dilution
Replicate 1	21.6	817.8	105.2	28.9	1529.8	366.7	<10	10 (50 for tartrate)
Replicate 2	21.5	817.3	105.6	28.7	1544.0	366.7	<10	10 (50 for tartrate)
Replicate 3	21.5	818.9	105.3	28.7	1528.8	366.7	<10	10 (50 for tartrate)
Average	21.5	818.0	105.3	28.8	1534.2	366.7	<10	
RSD (%)	0.040	0.098	0.19	0.44	0.56	0.009		

Chromatography


Calibration Standards, all analytes calibrated 1 – 100 ppm



Sulfite working standards prepared in 2% isopropanol




Red wine, 1:10 dilution, overlay of three injections



Red Wine (ppm)	Chloride	Phosphate	Malate	Sulfite	Tartrate	Sulfate	Oxalate
Average	60.0	771.6	92.1	26.7	1755.6	553.0	<10
RSD (%)	0.027	0.23	0.064	1.9	0.092	0.007	
Acetate not quantified.							

White wine, 1:10 dilution, overlay of three injections



White Wine (ppm)	Chloride	Phosphate	Malate	Sulfite	Tartrate	Sulfate	Oxalate
Average	21.5	818.0	105.3	28.8	1534.2	366.7	<10
RSD (%)	0.040	0.098	0.19	0.44	0.56	0.009	

Acetate not quantified.